

Helium pycnometry reference sheet


Phil McCausland, Earth Sciences

pmccausl@uwo.ca :: (519) 661-2111 x-88008



Instrument: Quantachrome Multipycnometer

$$V_{\text{fragment}} = V_{\text{sample cell}} - V_{\text{tank}} \left(\frac{P_{\text{initial}}}{P_{\text{final}}} - 1 \right)$$



Sample cells

Large volume	148.483 cm ³ ; depth 7.3 cm; inner diam. 4.85 cm
Medium volume	58.003 cm ³ ; depth 3.8 cm; inner diam. 3.95 cm
Small volume	28.827 cm ³ ; depth 3.6 cm; inner diam. 2.45 cm
Micro volume	12.364 cm ³ ; depth 2.3 cm; inner diam. 1.55 cm

Reference cells ("tank")

Large tank	88.595 cm ³
Small tank	11.629 cm ³
Micro tank	6.177 cm ³

Aluminum inserts (approx. volume)

Large A	51.491 cm ³
Large B	52.454 cm ³
Medium	40.367 cm ³
Small	11.281 cm ³
Micro	2.642 cm ³

Note: Aluminum inserts should be *measured* to obtain correct grain volume to subtract from Sample cell volume. Approx. volumes are from mass/Al density.

Reference spheres (NIST)

Large sphere	56.5592 ± 0.0023 cm ³
Medium sphere	28.9583 ± 0.0018 cm ³
Small sphere	7.0699 ± 0.0006 cm ³
Micro sphere	1.0725 ± 0.0004 cm ³